



PHILIPS

Computed tomography

Iconic **color** quantification

Philips IQon Spectral CT

At long last, color

The Philips IQon Spectral CT is the world's first and only spectral-detector CT, built from the ground up for spectral imaging. It delivers on-demand color quantification and the ability to characterize structures that's radically simple and low dose. Now, with the IQon Spectral CT, every scan can be spectral on demand.



Iconic quantification

Color quantification adds spectral resolution to your image quality delivering not just anatomical information but the ability to identify and characterize structures based on material content.

On-demand spectral analysis


Prospective and retrospective spectral results in one scan – without the need for special modes.

Iconic innovation

World's first spectral-detector CT built from the ground up for spectral imaging. And you get all this without complexity and at low dose. This is truly an iconic moment in the history of imaging. Your CT world is now in living color.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

Color **informs**. Color **quantifies**. Color **clarifies**.



Identify and quantify


Featuring first-of-its-kind technology, the IQon Spectral CT allows you to use color within CT images to identify the composition of what you see. Through this quantitative approach, you add spectral resolution to your image quality. So you not only get the anatomical information that you are used to with CT, but also uncover the ability to characterize tissues and structures based on material content. This is designed to help you improve your patient's care.



CT Angio Neck demonstrating different levels of MonoE.



Chest PE study evaluated with iodine overlay.



Chest PE study evaluated with Effective Z overlay.

Spectral on-demand made **simple**

Philips IQon Spectral CT allows for on-demand retrospective data analysis with your traditional workflow and without a special acquisition mode. This takes the guesswork out of multi-energy acquisitions, making it easy to use and allowing for routine spectral use. Retrospective spectral analysis is made possible through the iPatient platform, so you can experience spectral CT without the need for any special protocols.

Keep your traditional workflow

We designed Philips advances like iPatient with IQon Spectral CT in mind. With IQon Spectral CT, integrate the personalized quantitation of IQon Spectral CT directly into your scanning and reading workflow with full dose management.


You scan as you normally do and the spectral information is there, at your fingertips, when you need it. Now with the Philips IQon Spectral CT, every scan can be spectral on demand.




The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

Spectral is **always** on

With IQon Spectral CT, no prospective decision is required.
Enjoy spectral results anytime with the Spectral CT Viewer.



Advanced imaging tools allow for retrospective analysis that automatically generates spectral results with interactive analysis.



Results are ready to read on the PACS.


The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

Iconic innovation

With the launch of the IQon Spectral CT, we enrich the realm of clinical information and enable the “**and**” in CT. Through the uniqueness of the Philips detector-based spectral approach – **and** the NanoPanel Prism design – high and low energies live in the same time and space.


Advances in material science

The uniqueness of the Philips detector-based spectral approach and the NanoPanel Prism design allows you to get the conventional anatomical information that you are used to from your CT **and**, at the same time, get color quantification **and** the ability to characterize structures **and** monoenergetic image information. All in one scan, simultaneously. **And** you get it without increased complexity **and** at low dose.



NanoPanel Prism

- Low-dose, simultaneous spectral energy separation
- Low noise with Elite electronics
- No sensitivity to afterglow and no dead time



Yttrium-based scintillator

- Optimized for energy separation and with low image noise
- High light output at low energy
- Simultaneous detection in both time and space with negligible intra-layer scatter

Low dose spectral

Superb in black and white and now color, Philips is taking CT image quality at low dose to new heights with the introduction of IQon Spectral CT. Now you add spectral resolution to your image quality and maintain the high-quality images you have come to expect in noise reduction, spatial resolution, and low-contrast resolution.

Designed with Philips Spectral CT in mind




IMR* achieves 60–80% dose reduction and improves conventional image quality. These benefits are achieved simultaneously.

- 60–80% lower dose and IQ improvement
- 73–90% image noise reduction
- 2.5x–3.6x low contrast detectability
- 1.2x–1.7x spatial resolution



IMR off



IMR on

* In clinical practice, the use of IMR may reduce CT patient dose depending on the clinical task, patient size, anatomical location, and clinical practice. A consultation with a radiologist and a physicist should be made to determine the appropriate dose to obtain diagnostic image quality for the particular clinical task. Lower image noise, improved spatial resolution, improved low-contrast detectability, and/or dose reduction, were tested using reference body protocols. All metrics were tested on phantoms. Dose reduction assessments were performed using 0.8 mm slices, and tested on the MITA CT IQ Phantom (CCT183, The Phantom Laboratory), using human observers. Data on file.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

IQon Spectral CT image gallery

**Beam hardening
improvement in
the frontal brain**

Scan parameters

120 kVp


330 mAs

iDose⁴ – Level 3

Coverage – 16.8 cm

Scan time – 6.2 s

CTDI_{vol} – 58.8 mGy



Images courtesy of Hadassah Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.



120 kVp



80 keV

**Beam hardening demonstrated
on a MonoE 80 keV compared
to the conventional CT
of the sagittal brain**

Scan parameters


120 kVp

330 mAs

Coverage – 16.8 cm

Scan time – 10.8 s

CTDI_{vol} – 58.8 mGy



Images courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

CTA carotid with Spectral Results demonstrating the effects of different levels of MonoE on plaque

Scan parameters

120 kVp


250 mAs

iDose⁴ – Level 3

Coverage – 30.6 cm

Scan time – 3.3 s


CTDI_{vol} – 23.8 mGy



MonoE 50 keV



MonoE 75 keV



MonoE 100 keV

Images courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

**Step & Shoot Cardiac with
spectral results comparing
50 keV, 70 keV, and 90 keV**

Scan parameters

120 kVp

300 mAs

iDose⁴ – Level 3


Coverage – 16.1 cm

Scan time – 0.3 s


CTDI_{vol} – 37.5 mGy



50 keV



70 keV



90 keV

Images courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

Chest PE study evaluated with iodine overlay

Scan parameters

120 kVp

87 mAs

Coverage – 27.6 cm

Scan time – 5.6 s

CTDI_{vol} – 8.4 mGy

DLP – 284.8 mGy*cm



120 kVp and 45 keV chest with contrast comparison

Scan parameters

120 kVp

111 mAs

iDose⁴ – Level 3

Coverage – 29.0 cm

Scan time – 2.6 s

CTDI_{vol} – 10.5 mGy



Images courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

CTA comparing 120 kVp to 50 keV of the thoracic abdominal aorta

Scan parameters

120 kVp

130 mAs

iDose⁴ – Level 3

Coverage – 32.1 cm

Scan time – 0.21 s

CTDI_{vol} – 19.4 mGy




Image courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

**Routine CT abdomen
with spectral results
of liver metastasis**

Scan parameters


120 kVp

465 mAs

Coverage – 59.5 cm

Scan time – 26.9 s

CTDI_{vol} – 44.2 mGy



Images courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

CTA abdomen MonoE 70 keV demonstrating different spectral capabilities

Scan parameters

120 kVp




115 mAs

iDose⁴ – Level 3

Coverage – 43.4 cm

Scan time – 13.8 s

CTDI_{vol} – 10.9 mGy



Images courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

**IQon Spectral CT
abdomen and pelvis
images of IMR off
and IMR on**

Scan parameters

120 kVp

170 mAs

Coverage – 137.0 cm

Scan time – 25.6 s

CTDI_{vol} – 19.1 mGy



IMR off



IMR on



IMR off



IMR on

Images courtesy of Hadassah Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

**120 kVp and 50 keV
70 second injection delay**

Scan parameters

120 kVp



181 mAs

iDose⁴ – Level 4

Coverage – 69.0 cm

Scan time – 10.2 s

CTDI_{vol} – 16.2 mGy



Images courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

**CTA abdomen with
70-second injection
delay at MonoE
40 keV**

Scan parameters

120 kVp

242 mAs

iDose⁴ – Level 2

Coverage – 65.4 cm

Scan time – 9.6 s

CTDI_{vol} – 23.0 mGy




Image courtesy of University Hospitals Case Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

**Ankle results comparing
conventional image to
170 keV MonoE image**

Scan parameters

120 kVp

150 mAs

iDose⁴ – Level 3

Coverage – 25.9 cm

Scan time – 13.9 s

CTDI_{vol} – 16.7 mGy



Conventional



MonoE 170 keV

**Right hip results
comparing
conventional image to
150 keV MonoE image**

Scan parameters

120 kVp

270 mAs

iDose⁴ – Level 3

Coverage – 25.2 cm

Scan time – 7.7 s

CTDI_{vol} – 16.3 mGy



Conventional



MonoE 150 keV

Images courtesy of Hadassah Medical Center.

The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

Enhancing the capabilities of your existing iCT and Ingenuity CT family scanners, the SmartPath upgrade offers easy access to knowledge-based iterative reconstruction.



Optimize your system's performance both now and in the future with regular and ongoing updates, including functionality improvements and remote technical support.



Enhance your equipment with regular technology upgrades, and take advantage of the newest features and capabilities.



Transform your investment at the end of your system's life by transitioning seamlessly to a next-generation solution or refurbished option.

The images and descriptions contained herein provide technical specifications and optional features which may not be included with the standard system configuration. Contact your local Philips Representative for complete specific system details.

Some or all of the products, features, and accessories shown or described herein may not be available in your market. Please contact your local Philips Representative for availability.

CT performance specifications represent typical measured values.
The IQon Spectral CT is not yet CE Marked. Not available for sale in all regions.

© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.



Please visit www.philips.com/IQonSpectralCT

Printed in The Netherlands.
4522 991 05871 * DEC 2014